

Inactivation Type Virus Sampling Tube (Special For COVID Test)

Expected Usage :

For collection, transportation and storage of clinical samples

Package Spec :

50 tests/kit

Product Formation:

One OP swab + one sampling tube with 3 mL preservation liquid



Safe

Effectively inactivate virus to ensure biosafety



Stable

Special ingredients to ensure its integrity



Convenient

Allow room temperature storage of RNA

Operation Procedure

1



Use the swab to wipe the tonsils and the posterior pharyngeal wall, collect the sample.

2



Open the tube and put the swab into the tube.

3



Break off the swab and close the lid .

4



Record the sample information .